

| | |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name: | ABIS-1 |
| Description: | <p>CryoSafe ABIS-1 Class Liquid Nitrogen Tank Sentry Liquid Nitrogen Controller and Alarm features:</p> <ul style="list-style-type: none"> • Automatic LN2 monitoring and filling • Diagnostics at start-up • 4 thermistor array • LN2 level readout in LOW, NORMAL, and HIGH (direct sensor LN2 contacts) • Digital temperature display • Temp display units: C • Audible and visual alarm <ul style="list-style-type: none"> ○ LN2 level alarm (low and high level) ○ Sensor failure alarm ○ LN2 supply alarm ○ Power failure alarm (remote only) • Remote contacts • Controller surge protection • Low cost replacement <p>Model specifications: 10,400 (2ml) total vial capacity 235 total bag Capacity (250 ml Bags) 165 Liter liquid nitrogen capacity Static Evaporation Rate (Liters/Day): 5 Static Holding Time (Days): 33 Internal diameter (in/mm): 21/533.4 Overall height (in/mm): 44/1117.6 Usable height (in/mm): 29/736.6 Exterior Width (in/mm): 23.1/586.7 Exterior Height (in/mm): 30.5/774.7 Weight full (lbs/kg): 537/243.6 Shipping weight (lbs/kg): 400/181.4 Weight empty (lbs/kg): 245/111.1 Power Requirement: 24 VAC Three year parts and labor warranty Five year warranty on vacuum</p> |