

Product Name:	ABIS-2
Description:	<p>CryoSafe ABIS-2 Class Liquid Nitrogen Tank Sentry Liquid Nitrogen Controller and Alarm features:</p> <ul style="list-style-type: none"> • Automatic LN2 monitoring and filling • Diagnostics at start-up • 4 thermistor array • LN2 level readout in LOW, NORMAL, and HIGH (direct sensor LN2 contact) • Digital temperature display • Temp display units: C • Audible and visual alarm <ul style="list-style-type: none"> ○ LN2 level alarm (low and high level) ○ Sensor failure alarm ○ LN2 supply alarm ○ Power failure alarm (remote only) • Remote contacts • Controller surge protection • Low cost replacement <p>Model specifications: 24,050 (2ml) total vial capacity 535 total bag capacity (based on Fenwal #4R9953 Liquid nitrogen capacity (Liters): 356 Static Evaporation Rate (Liters/Day): 7 Static Holding Time (Days): 52 Internal diameter (in/mm): 31/787.4 Overall height (in/mm): 44/1117.6 Usable height (in/mm): 29/736.6 Exterior width (in/mm): 34/863.6 Exterior depth (in/mm): 38.5/977.9 Weight full (lbs/kg): 1046/474.5 Weight empty (lbs/kg): 405/183.7 Shipping weight (lbs/kg): 640/290.3 Power requirement: 24 VAC Three year parts and labor warranty Five year warranty on vacuum</p>