

| | |
|---------------|--|
| Product Name: | CM 3 PS |
| Description: | <p>CryoMizer Class Cryogenic Freezer</p> <p>CryoMizer Plastic Package System Includes:</p> <p>(5) Stainless Steel Racks, Including Center Max Rack</p> <p>(50) 81 or 100 Cell Plastic Boxes and Grids (Please specify when ordering)</p> <p>Rollerbase for CM-3</p> <p>One Pair CryoGloves</p> <p>Liquid Level Measuring Ruler</p> <p>Model Specifications:</p> <p>5,000 (2ml) Total Vial Capacity</p> <p>50 Box Capacity (100 or 81 Vials per box) (Boxes not included)</p> <p>5 Rack Capacity (10 boxes per rack)</p> <p>121 Liter Liquid Nitrogen Capacity</p> <p>Static Evaporation Rate (Liters/Days): 0.99</p> <p>Static Holding Time: 70 days</p> <p>Neck Diameter (in/mm): 8.5/215.9</p> <p>Exterior Diameter (in/mm): 22/558.8</p> <p>Exterior Height (in/mm): 38/965.2</p> <p>Weight Empty (lbs/kg): 81/36.7</p> <p>Weight Full (lbs/kg): 300/136.1</p> <p>Shipping Weight (lbs/kg) includes rollerbase: 210/95.3</p> <p>Three Year Parts and Labor Warranty</p> <p>Five Years Vacuum Warranty</p> |