

Product Name:	CM 3
Description:	<p>CryoMizer Class Cryogenic Freezer</p> <p>CryoMizer System Includes:</p> <p>(5) Stainless Steel Racks, Including Center Max Rack</p> <p>Rollerbase for CM-3</p> <p>Liquid Level Measuring Ruler</p> <p>Model Specifications:</p> <p>5,000 (2ml) Total Vial Capacity</p> <p>50 Box Capacity (100 or 81 Vials per box) (Boxes not included)</p> <p>5 Rack Capacity (10 boxes per rack)</p> <p>121 Liter Liquid Nitrogen Capacity</p> <p>Static Evaporation Rate (Liters/Days): 0.99</p> <p>Static Holding Time: 70 days</p> <p>Neck Diameter (in/mm): 8.5/215.9</p> <p>Exterior Diameter (in/mm): 22/558.8</p> <p>Exterior Height (in/mm): 38/965.2</p> <p>Weight Empty (lbs/kg): 81/36.7</p> <p>Weight Full (lbs/kg): 300/136.1</p> <p>Shipping Weight (lbs/kg) includes rollerbase: 170/77.1</p> <p>Three Year Parts and Labor Warranty</p> <p>Five Years Vacuum Warranty</p>