

Product Name:	SSC 21
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer</p> <p>Sample Storage in Canisters includes:</p> <p>9 canisters</p> <p>½ cc Straw (1 level bulk): 3870</p> <p>Liquid Nitrogen Capacity (Liters): 21</p> <p>Static Evaporation Rate (Liters/Day): 0.35</p> <p>Normal working time (days): 38</p> <p>Neck diameter (in/mm): 3.5/88.9</p> <p>Exterior height (in/mm): 17.2/436.8</p> <p>Exterior diameter (in/mm): 18.2/462.3</p> <p>Cannister Height (in/mm): 5/127</p> <p>Cannister Diameter (in/mm): 2.59/65.8</p> <p>Weight full (lbs/kg): 61.4/27.9</p> <p>Weight empty (lbs/kg): 24/10.9</p> <p>Includes liquid level measuring rod at NO CHARGE</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>