

Product Name:	SSC 47-6
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer</p> <p>Sample Storage in Canisters includes:</p> <p>6 canisters</p> <p>Total Vial Capacity (5/cane): 1320</p> <p>½ cc Straw (1 level bulk): 6200</p> <p>½ cc Straw (10 Cane): 3900+</p> <p>Liquid Nitrogen Capacity (Liters): 47.4</p> <p>Static Evaporation Rate (Liters/Day): 0.39</p> <p>Normal working time (days): 76</p> <p>Neck diameter (in/mm): 5/127</p> <p>Exterior height (in/mm): 26.5/673</p> <p>Exterior diameter (in/mm): 20/508</p> <p>Cannister Height (in/mm): 11/279</p> <p>Cannister Diameter (in/mm): 4/101.6</p> <p>Weight full (lbs/kg): 120.4/54.6</p> <p>Weight empty (lbs/kg): 42/19.1</p> <p>Includes liquid level measuring rod at NO CHARGE</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>