

Product Name:	SSC ET 33
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer with Extended Time Sample Storage in Canisters includes:</p> <p>6 canisters</p> <p>Total Vial Capacity (5/cane): 150</p> <p>½ cc straw capacity (1 level bulk): 780</p> <p>½ cc Straws (10 Cane): 540</p> <p>Liquid Nitrogen capacity (Liters): 33</p> <p>Static Evaporation Rate (Liters/Day): 0.13</p> <p>Normal working time (Days): 182</p> <p>Neck diameter (in/mm): 2.18/55.4</p> <p>Exterior height (in/mm): 25.9/657.9</p> <p>Exterior diameter (in/mm): 18.2/462.3</p> <p>Cannister Height (in/mm): 11/279.4</p> <p>Cannister Diameter (in/mm): 1.5/38.1</p> <p>Weight empty (lbs/kg): 34/15.4</p> <p>Weight full (lbs/kg): 93.4/42.4</p> <p>Includes liquid level measuring rod at NO CHARGE</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>