

Product Name:	VS 2
Description:	<p>CryoSafe Cryogenic Vapor Shipper Vapor Shipper includes: 1 canister ½ cc Straw (1 level bulk): 88 ½ cc Straws (10 Cane): 60 Total Vial Capacity (6/cane): 24 Liquid Nitrogen Capacity (Liters): 2.9 Static Evaporation Rate (Liters/Day): 0.20 Static Holding Time (Days): 14 Neck diameter (in/mm): 1.4/35.6 Exterior height (in/mm): 19.5/495.3 Exterior diameter (in/mm): 7.2/182.9 Cannister Height (in/mm): 11/279.4 Cannister Diameter (in/mm): 1.2/30.5 Weight full (lbs/kg): 15.5/7 Weight Empty (lbs/kg): 9/4.1 Weight Charged (lbs/kg): 13.7/6.2 Three year parts and labor warranty Three year vacuum warranty PSC 1 –Protective shipping carton available</p>