

Product Name:	VS 4
Description:	<p>CryoSafe Cryogenic Vapor Shipper Vapor Shipper includes: 7 canisters ½ cc Straw (1 level bulk): 4354 ½ cc Straws (10 Cane): 3080 Total Vial Capacity (6/cane): 1134 Liquid Nitrogen Capacity (Liters): 4.2 Static Evaporation Rate (Liters/Day): 0.35 Static Holding Time (Days): 12 Neck diameter (in/mm): 3.8/96.5 Exterior height (in/mm): 22/558.8 Exterior diameter (in/mm): 18.3/464.8 Cannister Height (in/mm): 11/279.4 Cannister Diameter (in/mm): 3.1/78.7 Weight full (lbs/kg): 87/39.5 Weight Empty (lbs/kg): 30.5/13.8 Weight Charged (lbs/kg): 40.9/18.6 Three year parts and labor warranty Three year vacuum warranty PSC 4 –Protective shipping carton available</p>