

Product Name:	VS 5
Description:	<p>CryoSafe Cryogenic Vapor Shipper  Vapor Shipper includes:  1 canister  ½ cc Straw (1 level bulk): 210  ½ cc Straws (10 Cane): 120  Total Vial Capacity (6/cane): 48  Liquid Nitrogen Capacity (Liters): 4.3  Static Evaporation Rate (Liters/Day): 0.20  Static Holding Time (Days): 21  Neck diameter (in/mm): 2/50.8  Exterior height (in/mm): 19.4/492.8  Exterior diameter (in/mm): 8.7/221  Cannister Height (in/mm): 11/279.4  Cannister Diameter (in/mm): 1.81/46  Weight full (lbs/kg): 21/9.5  Weight Empty (lbs/kg): 11.6/5.3  Weight Charged (lbs/kg): 19.1/8.7  Three year parts and labor warranty  Three year vacuum warranty  PSC 1 –Protective shipping carton available</p>