

Product Name:	VS 7
Description:	<p>CryoSafe Cryogenic Vapor Shipper</p> <p>Vapor Shipper includes:</p> <p>1 Rack</p> <p>Total Vial Capacity (6/cane): 500</p> <p>Total Bag Capacity (250 ml): 10</p> <p>Liquid Nitrogen Capacity (Liters): 8.5</p> <p>Static Evaporation Rate (Liters/Day): 0.85</p> <p>Static Holding Time (Days): 10</p> <p>Neck diameter (in/mm): 8.5/215.9</p> <p>Exterior height (in/mm): 21.5/546.1</p> <p>Exterior diameter (in/mm): 14.5/368.3</p> <p>Weight full (lbs/kg): 54.5/24.7</p> <p>Weight Empty (lbs/kg): 26.1/11.8</p> <p>Weight Charged (lbs/kg): 39/17.7</p> <p>Three year parts and labor warranty</p> <p>Three year vacuum warranty</p> <p>Includes protective shipping carton at NO CHARGE</p>