

Product Name:	VS 8
Description:	<p>CryoSafe Cryogenic Vapor Shipper Vapor Shipper includes: Total Vial Capacity (6/cane): 966 Bulk Total Bag Capacity (250 ml): 10 Liquid Nitrogen Capacity (Liters): 10 Static Evaporation Rate (Liters/Day): 0.70 Static Holding Time (Days): 14 Neck diameter (in/mm): 8.5/215.9 Exterior height (in/mm): 23/584.2 Exterior diameter (in/mm): 15/381 Weight full (lbs/kg): 73/33.1 Weight Empty (lbs/kg): 33/15 Weight Charged (lbs/kg): 49/22.2 Three year parts and labor warranty Three year vacuum warranty Includes protective shipping carton at NO CHARGE</p>