

Product Name:	VS 9
Description:	<p>CryoSafe Cryogenic Vapor Shipper Vapor Shipper includes: 6 Cannisters ½ cc Straw (1 level bulk): 780 ½ cc Straws (10 Cane): 540 Total Vial Capacity (6/cane): 180 Liquid Nitrogen Capacity (Liters): 12.3 Static Evaporation Rate (Liters/Day): 0.09 Static Holding Time (Days): 60 Neck diameter (in/mm): 2/50.8 Exterior height (in/mm): 25.7/652.8 Exterior diameter (in/mm): 14.5/368.3 Cannister Height (in/mm): 11/279.4 Cannister Diameter (in/mm): 1.5/38.1 Weight full (lbs/kg): 56.5/25.6 Weight Empty (lbs/kg): 25.7/11.6 Weight Charged (lbs/kg): 44.5/20.1 Three year parts and labor warranty Three year vacuum warranty</p>