

Product Name:	ASP-3 PS
Description:	<p>CryoSafe ASP Series Liquid Nitrogen Tank ASP-3 PS Package System includes: Stainless steel tank with cabinet Recessed casters Keyed lid lock Leveling platform 6 ft. stainless steel armor plated transfer hose (1) pair of Cryo gloves (28) large stainless steel racks (6) small stainless steel racks (364) 81/100 cell plastic boxes (78) 25 cell plastic boxes Liquid level measuring ruler Sentry plus controller Sentry Plus Auto-Fill Controller features:</p> <ul style="list-style-type: none"> • Power-Guard battery back up • Freeze-Guard overflow protection • Quick-chill temperature recovery • Auto-defog visibility safety • Splash-Guard false alarm prevention • Help button, user assistance • System log, alarm log, and temperature log with total 500,000 events • Serial port, 2 way communications • 4-20mA port • Ready for icSmart remote monitoring system • Keyless entry to power and control settings • Audio and visual alarms: <ul style="list-style-type: none"> ○ High temperature ○ High/low LN2 level ○ Valve stuck open ○ Lid open ○ Power failure ○ Sensor failure ○ Low LN2 supply ○ Unauthorized access ○ Usage warning ○ Temp calibration ○ Low battery • LN2 level readout +/- 1" (direct sensor LN2 contact) • LN2 supply monitor • LN2 usage calculation • Thermocouple auto-calibration • Temperature control • Digital temperature display

	<ul style="list-style-type: none">• Temp display units: C, F, K, R• Automatic filling• Diagnostics at start-up• Maintenance reminder• 8 thermistor array <p>Model specifications: 38,350 (2ml) total vial capacity 870 total bag capacity 590 liter liquid nitrogen capacity Static Evaporation Rate (liters/Day): 8 Static Holding Time: 74 Days Power Requirement: 16.5 VAC Exterior Diameter (in/mm): 42/1066.8 Overall Height (in/mm): 49/1244.6 Usable Height (in/mm): 29/736.6 Internal Diameter (in/mm): 39/990.6 Weight empty (lbs./kg): 565/256.3 Weight full (lbs./kg): 1616/733 Shipping weight (lbs./kg): 705/319.8 Three year parts and labor warranty Five year warranty on vacuum</p>
--	--