

Product Name:	CGSP-1
Description:	<p>The CryoGuard Class with Sentry Plus Auto-Fill Controller</p> <p>CryoGuard Stainless Steel Tank includes:</p> <p>Stainless steel outer &amp; inner tank wrapped with super insulate</p> <p>(5) lockable casters</p> <p>Hinged lid with locking provision</p> <p>Foldable skid free step</p> <p>Inner aluminum inventory carousel</p> <p>6 ft. armor plated transfer hose</p> <p>Fully operational battery back-up</p> <p>Sentry plus controller</p> <p>Sentry Plus Auto-Fill Controller features:</p> <ul style="list-style-type: none"> <li>• Power-Guard battery back up</li> <li>• Freeze-Guard overfill protection</li> <li>• Quick-chill temperature recovery</li> <li>• Auto-defog visibility safety</li> <li>• Splash-Guard false alarm prevention</li> <li>• Help button, user assistance</li> <li>• System log, alarm log, and temperature log with total 500,000 events</li> <li>• Serial port, 2 way communications</li> <li>• 4-20mA port</li> <li>• Ready for icSmart remote monitoring system</li> <li>• Keyless entry to power and control settings</li> <li>• Audio and visual alarms: <ul style="list-style-type: none"> <li>○ High temperature</li> <li>○ High/low LN2 level</li> <li>○ Valve stuck open</li> <li>○ Lid open</li> <li>○ Power failure</li> <li>○ Sensor failure</li> <li>○ Low LN2 supply</li> <li>○ Unauthorized access</li> <li>○ Usage warning</li> <li>○ Temp calibration</li> <li>○ Low battery</li> </ul> </li> <li>• LN2 level readout +/- 1" (direct sensor LN2 contact)</li> <li>• LN2 supply monitor</li> <li>• LN2 usage calculation</li> <li>• Thermocouple auto-calibration</li> <li>• Temperature control</li> <li>• Digital temperature display</li> <li>• Temp display units: C, F, K, R</li> <li>• Automatic filling</li> <li>• Diagnostics at start-up</li> <li>• Maintenance reminder</li> <li>• 8 thermistor array</li> </ul>

	<p>Model specifications:</p> <p>Total vial capacity (2ml): 19,500</p> <p>Total Bag Capacity (250 mL Bag): 500</p> <p>LN2 capacity (vapor phase): 36</p> <p>Total LN2 capacity (liters): 371</p> <p>Static evaporation rate (liters/day): 8</p> <p>Overall height (in/mm): 53/1346.2</p> <p>Reach over height (in/mm): 42.5/1079.5</p> <p>Usable height (in/mm): 30/762</p> <p>External diameter (in/mm): 34/863.6</p> <p>Internal diameter (in/mm): 29.5/749.3</p> <p>Neck Diameter (in/mm): 13/330.2</p> <p>Weight empty: 650 lbs</p> <p>Weight full: 1311 lbs.</p> <p>Shipping weight (lbs/kg): 1000/453.6</p> <p>Power requirements: 16.5 VAC</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>
--	--