

Product Name:	CGSP-2
Description:	<p>The CryoGuard Class with Sentry Plus Auto-Fill Controller CryoGuard Stainless Steel Tank includes: Stainless steel outer & inner tank wrapped with super insulate (5) lockable casters Hinged lid with locking provision Foldable skid free step Inner aluminum inventory carousel 6 ft. armor plated transfer hose Fully operational battery back-up Sentry plus controller Sentry Plus Auto-Fill Controller features:</p> <ul style="list-style-type: none"> • Power-Guard battery back up • Freeze-Guard overflow protection • Quick-chill temperature recovery • Auto-defog visibility safety • Splash-Guard false alarm prevention • Help button, user assistance • System log, alarm log, and temperature log with total 500,000 events • Serial port, 2 way communications • 4-20mA port • Ready for icSmart remote monitoring system • Keyless entry to power and control settings • Audio and visual alarms: <ul style="list-style-type: none"> ○ High temperature ○ High/low LN2 level ○ Valve stuck open ○ Lid open ○ Power failure ○ Sensor failure ○ Low LN2 supply ○ Unauthorized access ○ Usage warning ○ Temp calibration ○ Low battery • LN2 level readout +/- 1" (direct sensor LN2 contact) • LN2 supply monitor • LN2 usage calculation • Thermocouple auto-calibration • Temperature control • Digital temperature display • Temp display units: C, F, K, R • Automatic filling • Diagnostics at start-up • Maintenance reminder • 8 thermistor array

	<p>Model specifications: Total vial capacity (2ml): 41,600 Total Bag Capacity (250 ml Bag): 990 LN2 capacity (liters / vapor phase): 80 Total LN2 capacity (liters): 698 Static evaporation rate (liters/day): 9 Overall height (in/mm): 53/1346.2 Reach over height (in/mm): 42.5/1079.5 Usable height (in/mm): 30/762 Exterior diameter (in/mm): 45/1143 Internal diameter (in/mm): 40.5/1028.7 Neck Diameter (in/mm): 18/457.2 Weight empty (lbs/kg): 920/417.3 Weight full (lbs/kg): 2164/981.6 Shipping weight (lbs/kg) includes all components: 1500/680.4 Power requirements: 16.5 VAC Three year parts and labor warranty Five year vacuum warranty</p>
--	---