

Product Name:	CGSP-3
Description:	<p>The CryoGuard Class with Sentry Plus Auto-Fill Controller CryoGuard Stainless Steel Tank includes: Stainless steel outer & inner tank wrapped with super insulate (5) lockable casters Hinged lid with locking provision Foldable skid free step Inner aluminum inventory carousel 6 ft. armor plated transfer hose Fully operational battery back-up Sentry plus controller Sentry Plus Auto-Fill Controller features:</p> <ul style="list-style-type: none"> • Power-Guard battery back up • Freeze-Guard overflow protection • Quick-chill temperature recovery • Auto-defog visibility safety • Splash-Guard false alarm prevention • Help button, user assistance • System log, alarm log, and temperature log with total 500,000 events • Serial port, 2 way communications • 4-20mA port • Ready for icSmart remote monitoring system • Keyless entry to power and control settings • Audio and visual alarms: <ul style="list-style-type: none"> ○ High temperature ○ High/low LN2 level ○ Valve stuck open ○ Lid open ○ Power failure ○ Sensor failure ○ Low LN2 supply ○ Unauthorized access ○ Usage warning ○ Temp calibration ○ Low battery • LN2 level readout +/- 1" (direct sensor LN2 contact) • LN2 supply monitor • LN2 usage calculation • Thermocouple auto-calibration • Temperature control • Digital temperature display • Temp display units: C, F, K, R • Automatic filling • Diagnostics at start-up • Maintenance reminder

	<ul style="list-style-type: none">• 8 thermistor array <p>Model specifications: Total vial capacity (2ml): 79,300 LN2 capacity (liters / vapor phase): 137 Total LN2 capacity (liters): 1269 Total bag capacity: 1920 Static evaporation rate (liters/day): 15 Overall height (in/mm): 53/1346.2 Reach over height (in/mm):42.5/1079.5 Usable height (in/mm): 30/762 External diameter (in/mm): 59.5/1508.8 Internal diameter (in/mm): 55/1397 Neck diameter: 24.5/622 Weight empty (lbs/kg): 1550/703.1 Weight full (lbs/kg) 3811/1728.2 Shipping weight (lbs/kg) includes all components: 2,000/907.2 Power requirements: 16.5 VAC Three year parts and labor warranty Five year vacuum warranty</p>
--	--