

Product Name:	CM 1 PS
Description:	<p>CryoMizer Class Cryogenic Freezer CryoMizer Plastic Package System Includes: (7) Stainless Steel Racks, including center Max Rack (35) 25 Cell Plastic Boxes and Grids Rollerbase for CM-1 One Pair CryoGloves Liquid Level Measuring Ruler Model Specifications: 875 (2ml) Total Vial Capacity 35 Box Capacity (25 vials per box) (boxes not included) 7 Rack Capacity (5 boxes per rack) 47.4 Liter Liquid Nitrogen Capacity Static Evaporation Rate (Liters/Day): 0.39 Static Holding time: 76 days Neck Diameter (in/mm): 5/127 Exterior Height (in/mm): 26.5/673.1 Exterior Diameter (in/mm): 20/508 Weight Empty (lbs/kg): 42/19.1 Weight Full (lbs/kg): 120.4/54.6 Shipping Weight (lbs/kg) - includes rollerbase: 85/38.5 Three Year Parts and Labor Warranty Five Years Vacuum Warranty</p>