

Product Name:	CM 2 CS
Description:	<p>CryoMizer Class Cryogenic Freezer</p> <p>CryoMizer Cardboard Package System Includes:</p> <p>(5) Stainless Steel Racks, Including Center Max Rack</p> <p>(25) 100 or 81 Cell Cardboard Boxes and Grids (Please specify when ordering)</p> <p>Rollerbase for CM-2</p> <p>One Pair CryoGloves</p> <p>Liquid Level Measuring Ruler</p> <p>Model Specifications:</p> <p>2,500 (2ml) Total Vial Capacity</p> <p>25 Box Capacity (100 or 81 Vials per box) (Boxes not included)</p> <p>5 Rack Capacity (5 boxes per rack)</p> <p>61 Liter Liquid Nitrogen Capacity</p> <p>Static Evaporation Rate (Liters/Days): 0.85</p> <p>Static Holding Time: 38 days</p> <p>Inner Neck Diameter (in/mm): 8.5"/215.9</p> <p>Exterior Diameter (in/mm): 22"/558.8</p> <p>Exterior Height (in/mm): 27.25"/692.2</p> <p>Weight Empty (lbs/kg): 58/26.3</p> <p>Weight Full (lbs/kg): 135/61.2</p> <p>Shipping Weight (lbs/kg) – includes rollerbase: 155/70.3</p> <p>Three Year Parts and Labor Warranty</p> <p>Five Year Vacuum Warranty</p>