

Product Name:	CM 4 CS
Description:	<p>CryoMizer Class Cryogenic Freezer</p> <p>CryoMizer Cardboard Package System Includes:</p> <p>(7) Stainless Steel Racks, Including Center Max Rack</p> <p>(70) 100 or 81 Cell Cardboard Boxes and Grids (Please specify)</p> <p>Rollerbase for CM-4</p> <p>One Pair CyroGloves</p> <p>Liquid Level Measuring Ruler</p> <p>Model Specifications:</p> <p>7,000 Total (2ml) Vial Capacity</p> <p>70 Box Capacity (100 or 81 vials per box) (Boxes not included)</p> <p>7 Rack Capacity (10 boxes per rack)</p> <p>LN2 Capacity (Liters): 175</p> <p>Static Evaporation Rate (Liters/Day): 0.99</p> <p>Static Working Time: 104 Days</p> <p>Neck Diameter (in/mm): 8.5/215.9</p> <p>Exterior Diameter (in/mm): 26/660.4</p> <p>Exterior Height (in/mm): 37.75/958.8</p> <p>Weight Empty (lbs/kg): 103/46.7</p> <p>Weight Full (lbs/kg): 425/192.7</p> <p>Shipping Weight (lbs/kg) includes rollerbase: 210/95.3</p> <p>Three Year Parts and Labor Warranty</p> <p>Five Year Vacuum Warranty</p>