

Product Name:	LD 10
Description:	<p>CryoSafe Liquid Nitrogen Dewar</p> <p>Model specifications:</p> <p>Liquid nitrogen capacity (Liters): 10</p> <p>Static evaporation rate (Liters/Day): 0.18</p> <p>Neck diameter (in/mm): 2.18/55.4</p> <p>Exterior height (in/mm): 21.5/546.1</p> <p>Exterior diameter (in/mm): 10.3/261.6</p> <p>Weight full (lbs/kg): 31/14.1</p> <p>Weight empty (lbs/kg): 12/5.4</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>