

Product Name:	LD 20
Description:	<p>CryoSafe Liquid Nitrogen Dewar</p> <p>Model specifications:</p> <p>Liquid nitrogen capacity (Liters): 20</p> <p>Static evaporation rate (Liters/Day): 0.18</p> <p>Neck diameter (in/mm): 2.2/56</p> <p>Exterior height (in/mm): 24.5/622.3</p> <p>Exterior diameter (in/mm): 14.5/368.3</p> <p>Weight full (lbs/kg): 55/24.9</p> <p>Weight empty (lbs/kg): 19/8.6</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>