

Product Name:	LD 4
Description:	CryoSafe Liquid Nitrogen Dewar Model specifications: Liquid nitrogen capacity (Liters): 4 Static evaporation rate (Liters/Day): 0.19 Neck diameter (in/mm): 1.4/35 Exterior height (in/mm): 16.8/426.7 Exterior diameter (in/mm) 7.3/185.5 Weight full (lbs/kg): 13/5.9 Weight empty (lbs/kg): 6/2.7 Three year parts and labor warranty Five year vacuum warranty