

Product Name:	LD 50
Description:	<p>CryoSafe Liquid Nitrogen Dewar</p> <p>Model specifications:</p> <p>Liquid nitrogen capacity (Liters): 50</p> <p>Static evaporation rate (Liters/Day): 0.49</p> <p>Neck diameter (in/mm): 2.5/63.5</p> <p>Exterior height (in/mm): 30.5/774.7</p> <p>Exterior diameter (in/mm): 17/431.8</p> <p>Weight full (lbs/kg): 120/54.4</p> <p>Weight empty (lbs/kg): 31/14.1</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>