

Product Name:	SSC 20
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer</p> <p>Sample Storage in Canisters includes:</p> <p>6 canisters</p> <p>Total Vial Capacity (5/cane): 210</p> <p>½ cc Straw (1 level bulk): 1122</p> <p>½ cc Straws (10 Cane): 720</p> <p>Liquid Nitrogen Capacity (Liters): 20.5</p> <p>Static Evaporation Rate (Liters/Day): 0.09</p> <p>Normal working time (days): 135</p> <p>Neck diameter (in/mm): 2.18/55.4</p> <p>Exterior height (in/mm): 26/660</p> <p>Exterior diameter (in/mm): 14.5/368.3</p> <p>Cannister Height (in/mm): 11/279.4</p> <p>Cannister Diameter (in/mm): 1.65/41.9</p> <p>Weight full (lbs/kg): 58.4/26.5</p> <p>Weight empty (lbs/kg): 22/10</p> <p>Includes liquid level measuring rod at NO CHARGE</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>