

Product Name:	SSC 32
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer</p> <p>Sample Storage in Canisters includes:</p> <p>9 canisters</p> <p>Total Vial Capacity (5/cane): 855</p> <p>½ cc Straw (1 level bulk): 3960</p> <p>½ cc Straws (10 Cane): 2520</p> <p>Liquid Nitrogen Capacity (Liters): 32</p> <p>Static Evaporation Rate (Liters/Day): 0.35</p> <p>Normal working time (days): 57</p> <p>Neck diameter (in/mm): 3.81/96.7</p> <p>Exterior height (in/mm): 21.5/546.1</p> <p>Exterior diameter (in/mm): 18.2/462.3</p> <p>Cannister Height (in/mm): 11/279.4</p> <p>Cannister Diameter (in/mm): 2.62/66.5</p> <p>Weight full (lbs/kg): 87/39.5</p> <p>Weight empty (lbs/kg): 30/13.6</p> <p>Includes liquid level measuring rod at NO CHARGE</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>