

Product Name:	SSC 33
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer</p> <p>Sample Storage in Canisters includes:</p> <p>6 canisters</p> <p>Total Vial Capacity (5/cane): 360</p> <p>½ cc Straw (1 level bulk): 1764</p> <p>½ cc Straws (10 Cane): 1260</p> <p>Liquid Nitrogen Capacity (Liters): 33.4</p> <p>Static Evaporation Rate (Liters/Day): 0.14</p> <p>Normal working time (days): 154</p> <p>Neck diameter (in/mm): 2.75/69.8</p> <p>Exterior height (in/mm): 26/660.4</p> <p>Exterior diameter (in/mm): 18.2/464.3</p> <p>Cannister Height (in/mm): 11/279.4</p> <p>Cannister Diameter (in/mm): 2.22/56.4</p> <p>Weight full (lbs/kg): 94/42.6</p> <p>Weight empty (lbs/kg): 34/15.4</p> <p>Includes liquid level measuring rod at NO CHARGE</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>