

Product Name:	SSC ET 11
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer with Extended Time Sample Storage in Canisters includes:</p> <p>6 canisters</p> <p>Total Vial Capacity (5/cane): 150</p> <p>½ cc straw capacity (1 level bulk): 879</p> <p>½ cc Straws (10 Cane): 660</p> <p>Liquid Nitrogen capacity (Liters): 11</p> <p>Static Evaporation Rate (Liters/Day): 0.16</p> <p>Normal working time (Days): 43</p> <p>Neck diameter (in/mm): 2.18/55.4</p> <p>Exterior height (in/mm): 21.6/548.6</p> <p>Exterior diameter (in/mm): 10.2/259.1</p> <p>Cannister Height (in/mm): 11/279.4</p> <p>Cannister Diameter (in/mm): 1.65/41.9</p> <p>Weight empty (lbs/kg): 17/7.7</p> <p>Weight full (lbs/kg): 36.6/16.6</p> <p>Includes liquid level measuring rod at NO CHARGE</p> <p>Three year parts and labor warranty</p> <p>Five year vacuum warranty</p>