

Product Name:	SSC ET 3
Description:	<p>CryoSafe SSC Class Liquid Nitrogen Freezer with Extended Time Sample Storage in Canisters includes:</p> <ul style="list-style-type: none"> 6 canisters ½ cc straw capacity (1 level bulk): 1122 Liquid Nitrogen capacity (Liters): 3.6 Static Evaporation Rate (Liters/Day): 0.13 Normal working time (Days): 17 Neck diameter (in/mm): 2.18/55.4 Exterior height (in/mm): 16/406.4 Exterior diameter (in/mm): 8.7/220.9 Cannister Height (in/mm): 5/127 Cannister Diameter (in/mm): 1.65/41.9 Weight empty (lbs/kg): 8/3.6 Weight full (lbs/kg): 14.4/6.5 Includes liquid level measuring rod at NO CHARGE Three year parts and labor warranty Five year vacuum warranty