

Product Name:	VS 1
Description:	<p>CryoSafe Cryogenic Vapor Shipper Vapor Shipper includes: 1 canister ½ cc Straw (1 level bulk): 88 Liquid Nitrogen Capacity (Liters): 1.5 Static Evaporation Rate (Liters/Day): 0.19 Static Holding Time (Days): 8 Neck diameter (in/mm): 1.4/35.6 Exterior height (in/mm): 13.5/342.9 Exterior diameter (in/mm): 7.25/184.2 Cannister Height (in/mm): 5/127 Cannister Diameter (in/mm): 1.2/30.5 Weight full (lbs/kg): 9/4.1 Weight Empty (lbs/kg): 6/2.7 Weight Charged (lbs/kg): 8/3.6 Three year parts and labor warranty Three-year vacuum warranty PSC 1 –Protective shipping carton available</p>