

|               |   |
|---------------|---|
| Product Name: | VS 3  |
| Description:  | <p>CryoSafe Cryogenic Vapor Shipper<br/> Vapor Shipper includes:<br/> 1 canister<br/> ½ cc Straw (1 level bulk): 440<br/> ½ cc Straws (10 Cane): 280<br/> Total Vial Capacity (6/cane): 102<br/> Liquid Nitrogen Capacity (Liters): 3.6<br/> Static Evaporation Rate (Liters/Day): 0.26<br/> Static Holding Time (Days): 13<br/> Neck diameter (in/mm): 2.75/69.9<br/> Exterior height (in/mm): 18.4/467.4<br/> Exterior diameter (in/mm): 8.7/221<br/> Cannister Height (in/mm): 11/279.4<br/> Cannister Diameter (in/mm): 2.62/66.6<br/> Weight full (lbs/kg): 20/9.1<br/> Weight Empty (lbs/kg): 10/4.5<br/> Weight Charged (lbs/kg): 17.1/7.7<br/> Three year parts and labor warranty<br/> Three year vacuum warranty<br/> PSC 1 –Protective shipping carton available</p> |