

Product Name:	VS 6
Description:	<p>CryoSafe Cryogenic Vapor Shipper</p> <p>Vapor Shipper includes:</p> <p>5 canisters</p> <p>½ cc Straw (1 level bulk): 3750</p> <p>½ cc Straws (10 Cane): 2500</p> <p>Total Vial Capacity (6/cane): 840</p> <p>Liquid Nitrogen Capacity (Liters): 6.2</p> <p>Static Evaporation Rate (Liters/Day): 0.35</p> <p>Static Holding Time (Days): 16</p> <p>Neck diameter (in/mm): 3.81/96.7</p> <p>Exterior height (in/mm): 22/559</p> <p>Exterior diameter (in/mm): 14.5/368.3</p> <p>Cannister Height (in/mm): 11/279.4</p> <p>Cannister Diameter (in/mm): 3.2/81.3</p> <p>Weight full (lbs/kg): 60.2/27.3</p> <p>Weight Empty (lbs/kg): 25.1/11.3</p> <p>Weight Charged (lbs/kg): 36/16.3</p> <p>Three year parts and labor warranty</p> <p>Three year vacuum warranty</p> <p>PSC-3 – Protective Shipping Carton Available</p>